

ANOTHER SOLUTIONS

# APK 2500

## MULTI VALVE CONTROLLER

### MULTI VALVE CONTROLLER

APK2500 can control up to 4 valves independently. Gas or liquid sample can be available. This can be connected with all GC(Gas Chromatograph). Easy to setup, control via GC app or its own app(optional). This is useful to inject multi samples into one instrument.

### FEATURES

- Gas or liquid sample injection
- Minimized dead volume
- Heating Transfer line to prevent contamination
- With an instrument, fully automated sequence available
- Variable applications are ready via changing of valve paths

### SPECIFICATION

SPECIFICATIONS	
Temp. (2 Point)	Max. 200°C
Valve	VICI Valve (Micro Electrical or Air Actuator)
Dimension	200 x 500 x 530 (W x D x H mm)
Power	220 VAC, 50/60Hz, 0.7KW, 3A

